

AIR_EC_Capabilities for Milvue Suite

System name

OTHER_MILVUE_MVS

AI Algorithm -

1 I. SmartXpert

AI Algorithm Description

The AI algorithm is capable of analyzing X-ray images for performing automatic 2D measurements of orthopedic static parameters and bone age assessment. SmartXpert is indicated for the following anatomical regions and views:

- Left Hand: frontal view
- Shoulder: glenoid view also known as a 'true AP' or a 'Grashey view'
- Spine: full spine frontal and lateral views
- Pelvis: frontal view, weight-bearing
- Hip: false-profile view, weight-bearing
- Leg: Face and lateral Long Leg Radiographs, weight-bearing
- Foot: frontal and lateral views, weight-bearing

DICOM IODs implemented

- Comprehensive 3D SR Storage IOD
- Secondary Capture Image IOD

Result primitives encoded

- 6.5.3.2 Qualitative Findings
- 6.5.3.3 Measurements
- 6.5.3.4 Locations
- 6.5.3.7 Tracking Identifiers

Additional Information

Original Image



Results



AI Algorithm -

II. SmartUrgences

AI Algorithm Description

The AI algorithm is capable of analyzing X-ray images for detecting a variety of pathologies for Musculoskeletal (MSK) disorders and chest. For MSK, it identifies fractures in the following anatomical regions:

- Chest (including ribs),
- Clavicle,
- Shoulder,
- Arm, elbow, forearm, wrist, hand, fingers, hip, pelvis, femur, knee, leg (including tibia and fibula), ankle, foot

SmartUrgences is trained to detect the following anomalies:

- 1 Bone fracture
- 2 Joint dislocation
- 3 Elbow joint effusion
- 4 Pleural effusion
- 5 Pulmonary opacity
- 6 Pulmonary nodule
- 7 Pneumothorax

DICOM IODs implemented


- Comprehensive 3D SR Storage IOD
- Secondary Capture Image IOD

Result primitives encoded

- 6.5.3.2 Qualitative Findings
- 6.5.3.3 Measurements
- 6.5.3.4 Locations
- 6.5.3.7 Tracking Identifiers

Additional Information

Original Image



Results

